



READY ANKLE BY MOBILITY BRACES

Easy application ankle brace providing stability and support!

Are your patients **READY** to increase their mobility, comfort and confidence? The **READY ANKLE** is designed to work with the patient's shoe to help support and stabilize the ankle while assisting in providing confidence and pain relief for multiple indications and conditions. It is recommended for patients with various acute and chronic ankle injuries and/or conditions for use during activities of daily living.

FEATURES

- Simple 3-strap system makes application intuitive, quick and easy
- Lower circumferential strap wraps completely around the leg just above the ankle to provide improved compression and support
- Rigid medial and lateral upright struts help prevent abnormal inversion and eversion during daily activity
- Integrated footplate/functional hinge tracks with the user's ankle to allow for a more comfortable and natural movement
- Strap pull-tabs make application of the straps easier for patients with poor dexterity in their hands

INDICATIONS

Helps provide functional support and stability for:

- Chronic ankle instability
- Acute and chronic ankle sprains
- Prophylactic use
- Syndesmosis sprain
- Transition from walker boot or cast to daily activity
- Ankle injury rehabilitation of grade I, II and high ankle sprain and strains
- Provide additional inversion/eversion control for unstable ankles

ORDERING

RIGHT Part #	LEFT Part #	Size	Shoe Size (US)
MA01RSM	MA01LSM	Small	Men's size: up to 9 Women's size: up to 11
MA01RMD	MA01LMD	Medium	Men's size: 9 - 12 Women's size: up to 11 - 13.5
MA01RLG	MA01LLG	Large	Men's size: 12 and above Women's size: 13.5 and above

Sizing Notes: Sizing is based on US men's and women's shoe sizes.
Sizing may need to be adjusted for patients with disproportionate foot & ankle widths.

