



READY KNEE BY MOBILITY BRACES

Easy application knee brace providing stability and pain relief!

Are your patients **READY** to increase their mobility and decrease their pain? The **READY KNEE** is designed to help stabilize the knee and provide pain relief for various indications. Adjustable hinges allow for easy range-of-motion (ROM) control enabling the physician or practitioner the ability to set specific limits on flexion and extension. The **READY KNEE** comes with sixteen metal stops and includes a tool to make adjusting the ROM easy, providing maximum control and adjustability.

FEATURES

- Breathable fabric provides more airflow and increased comfort for better patient compliance
- Dual range-of-motion (ROM) hinges offer increased control and adjustability
- Easy-pull strap tabs make grabbing and pulling straps easier for patients with limited dexterity
- Wrap-around design with anterior closure allows for a customized fit and easier application than pull-on sleeves
- Rigid hinge arms can be bent to conform to patients anatomy providing a more customized and functional fit

INDICATIONS

Helps provide pain relief and functional support and stability for:

- Mild to moderate sprains and strains
- Medial/lateral ligament support
- Mild fractures
- Osteoarthritis
- Rheumatoid arthritis
- Meniscal cartilage derangement

ORDERING



Part Number	Size	Measurement
MK01UXS	XSmall	Circumference: < 15"
MK01USM	Small	Circumference: 15" to 17"
MK01UMD	Medium	Circumference: 17" to 20"
MK01ULG	Large	Circumference: 20" to 23.5"
MK01UXL	XLarge	Circumference: 23.5" to 26"
MK01U2X	2XLarge	Circumference: 26" to 29"
MK01U3X	3XLarge	Circumference: 29" to 34"
MK01U4X	4XLarge	Circumference: > 34"

Circumferential measurement of thigh taken 6" above knee center.