

# READY BACK BY MOBILITY BRACES

Easy application back brace providing support and pain relief!

Are your patients **READY** to increase their mobility and decrease their pain? The **READY BACK** is a **Lumbar Sacral Orthosis (LSO)** that provides direct support to the lower back and targets a wide variety of diagnoses. It produces intracavitary pressure to reduce load on intervertebral discs helping provide support from the sacrococcygeal junction to T-9 vertebrae, helping to relieve lower back pain and provide spinal support through comfortable immobilization. This is achieved with two removable lateral support panels along with an two removeable anterior (abdominal) panels and posterior (back) panel that work in conjunction with a patient controlled compression closure system. This system gently but firmly stabilizes the brace, allowing pressure to be relieved from associated discs or the injured area, resulting in a reduction in pain and improved posture.

## FEATURES

- Breathable fabric provides more airflow and increased comfort for better patient compliance
- Removable abdominal panels allow for increased or decreased levels of immobilization
- Zippered palm pockets on belt and thumb loops on closure system straps make application easier for patients with limited dexterity

## INDICATIONS

- Postural support
- Chronic Low Back Pain
- Herniated Discs
- Spinal Stenosis
- Spondylosis
- Spondylolisthesis
- Osteoarthritis
- Lumbago
- Muscle weakness

## ORDERING



Part Number	Size	Waist Circumference
MB01UXS	XSmall	Circumference: 23" to 26"
MB01USM	Small	Circumference: 26" to 30"
MB01UMD	Medium	Circumference: 30" to 35"
MB01ULG	Large	Circumference: 35" to 40"
MB01UXL	XLarge	Circumference: 40" to 44"
MB01U2X	2XLarge	Circumference: 44" to 48"
MB01U3X	3XLarge	Circumference: 48" to 53"
MB01U4X	4XLarge	Circumference: 53" to 59"
MB01U5X	5XLarge	Circumference: 59" to 65"