



READYlite BACK

BY MOBILITY BRACES

Lightweight, low-profile back brace providing support and pain relief!

Are your patients **READY** to increase their mobility and decrease their pain? The **READYlite BACK** is a **Lumbar Orthosis (LO)** that provides support from the L-1 to below the L-5 vertebra and targets a wide variety of diagnoses. It was designed to help relieve lower back pain and provide spinal support through comfortable immobilization. This is achieved with two removeable anterior (abdominal) panels and posterior (back) panel that work in conjunction with a patient controlled compression closure system. This system gently but firmly stabilizes the brace, allowing pressure to be relieved from associated discs or the injured area, resulting in a reduction in pain and improved posture.

FEATURES

- Breathable fabric provides more airflow and increased comfort for better patient compliance
- Zippered palm pockets on belt and thumb loops on closure system straps make application easier for patients with limited dexterity
- Patient controlled closure system makes adjusting the level of compression and support easy and quick

INDICATIONS

- Postural Support
- Chronic Low Back Pain
- Herniated Discs
- Spinal Stenosis
- Spondylosis
- Spondylolisthesis

ORDERING



Part Number	Size	Waist Circumference
MB02UXS	XSmall	Circumference: < 26"
MB02USM	Small	Circumference: 26" to 30"
MB02UMD	Medium	Circumference: 30" to 34"
MB02ULG	Large	Circumference: 34" to 38"
MB02UXL	XLarge	Circumference: 38" to 42"
MB02U2X	2XLarge	Circumference: 42" to 48"
MB02U3X	3XLarge	Circumference: 48" to 54"
MB02U4X	4XLarge	Circumference: > 54"