



READYlite KNEE

BY MOBILITY BRACES

Lightweight, low-profile knee brace providing stability and pain relief!

Are your patients **READY** to increase their mobility and decrease their pain? The **READYlite KNEE** is designed to help provide warmth, compression and stability of the knee while assisting in pain relief for various mild-to-moderate indications. The anterior closure design makes application quick and the monocentric hinges are lightweight and low-profile leading to better patient comfort and compliance.

FEATURES

- Monocentric hinges offer lightweight, low-profile medial/lateral stability and support
- Breathable fabric provides more airflow and increased comfort for better patient compliance
- Easy-pull strap tabs make grabbing and pulling straps easier for patients with limited dexterity
- Wrap-around design with anterior closure allows for a more customized fit and easier application than pull-on sleeves
- Open popliteal design prevents bunching behind the knee
- Hinges are removable for washing/maintenance

INDICATIONS

Helps provide pain relief and functional support and stability for:

- Mild to moderate sprains and strains of the knee
- Mild to moderate ligament instabilities of the knee
- Mild osteoarthritis of the knee

ORDERING



Part Number	Size	Measurement
MK03UXS	XSmall	Circumference: < 15"
MK03USM	Small	Circumference: 15" to 16.5"
MK03UMD	Medium	Circumference: 16.5" to 19.5"
MK03ULG	Large	Circumference: 19.5" to 23"
MK03UXL	XLarge	Circumference: 23" to 25.5"
MK03U2X	2XLarge	Circumference: 25.5" to 28"
MK03U3X	3XLarge	Circumference: 28" to 31"
MK03U4X	4XLarge	Circumference: > 31"

Circumferential measurement of thigh taken 6" above knee center.

