

KNEE BRACE SUSPENSION SLEEVE

Dual application suspension sleeve designed to help improve knee brace suspension, placement and comfort!

Are your patients having difficulty with their current knee brace? Are they experiencing fit and suspension issues? Are they dealing with comfort concerns? If they are, then that may lead to compliance concerns for you.

Are you **READY** to offer a solution to your patients for these issues? The Knee Suspension Sleeve is designed to aid in the proper placement, suspension and comfort of knee braces. Designed with dual application options, the Knee Suspension Sleeve can be applied either under or over the brace (based on patient preference) to improve fit, reduce brace migration and increase patient comfort.

FEATURES

- Breathable, perforated fabric provides more airflow and increased comfort for better patient compliance
- Dual application options provide 2 completely different ways to improve suspension and fit
- Easy-pull finger cut-outs make grabbing and pulling sleeve easier for patients with limited dexterity
- Trimmable base material allows for a more customized fit

INDICATIONS

Used in conjunction with a knee brace to provide improved suspension and comfort for patients needing stability for:

- Mild to moderate sprains and strains
- Medial/lateral ligament support
- Mild fractures
- Osteoarthritis
- Rheumatoid arthritis
- Meniscal cartilage derangement

ORDERING



| Part Number | Size | Knee Brace Sizes |
|-------------|------|--------------------|
| MK02UNI | OSFA | Works with XS - 4X |



OVER THE BRACE



UNDER THE BRACE