



READY KNEE BY MOBILITY BRACES

Easy application knee brace providing stability and pain relief!

Are your patients **READY** to increase their mobility and decrease their pain? The **READY KNEE** is designed to help stabilize the knee and provide pain relief for various indications. Adjustable hinges allow for easy range-of-motion (ROM) control enabling the physician or practitioner the ability to set specific limits on flexion and extension. The **READY KNEE** comes with sixteen metal stops and includes a tool to make adjusting the ROM easy, providing maximum control and adjustability.

FEATURES

- Breathable fabric provides more airflow and increased comfort for better patient compliance
- Dual range-of-motion (ROM) hinges offer increased control and adjustability
- Easy-pull strap tabs make grabbing and pulling straps easier for patients with limited dexterity
- Wrap-around design with anterior closure allows for a customized fit and easier application than pull-on sleeves
- Rigid hinge arms can be bent to conform to patients anatomy providing a more customized and functional fit

INDICATIONS

Helps provide pain relief and functional support and stability for:

- Mild to moderate sprains and strains
- Medial/lateral ligament support
- Mild fractures
- Osteoarthritis
- Rheumatoid arthritis
- Meniscal cartilage derangement

ORDERING



Part Number	Size	Measurement
MK01UXS	XSmall	Circumference: 14" to 15"
MK01USM	Small	Circumference: 15" to 16"
MK01UMD	Medium	Circumference: 16" to 18"
MK01ULG	Large	Circumference: 18" to 21"
MK01UXL	XLarge	Circumference: 21" to 24"
MK01U2X	XXLarge	Circumference: 24" to 27"
MK01U3X	XXXLarge	Circumference: 27" to 29"
MK01U4X	XXXXLarge	Circumference: 29" to 31"

Circumferential measurement of thigh 6" above knee.